

P1 Transduction

Prepare P1 Lysate

1. Grow o/n culture of strain you want to make P1 lysate of (has selective marker).
2. Dilute o/n culture 1/100 in 10 ml LB + 5 mM CaCl_2 + 0.2% glucose.
3. Grow 1 h at 37°C, no shaking. For strains that have to be grown at 32 ° C, grow for 1.5 h.
4. Add 200 μl of P1 lysate (IN THE HOOD) from stock**
5. Incubate at 37°C with shaking for 3 h or until lysis occurs (typically >2 h).
6. Transfer to 15 ml polypropylene tubes (IN THE HOOD) and add 200 μl CHCl_3 . Mix by vortexing. N.B. do not use polystyrene tubes as these will react with the CHCl_3 .
7. Centrifuge at full speed for 5 minutes.
8. Collect lysate and store with 200 μl CHCl_3 . Use in P1 transduction.

**To make more P1 lysate, amplify on strain MG1655 or another “wt” strain

P1 Transduction

- 1.** Label 4 eppendorf tubes # 1-4. Add 0.5 ml LB + 15mM CaCl₂ + 30mM MgCl₂ to each. Add 0.5 ml LB to tubes #2-4.
- 2.** Add 0.5 ml P1 lysate to tube #1. Transfer 100 µl from tube #1 to tube #2. Mix. Transfer 100 µl from tube #2 to tube #3. Mix.
- 3.** Add 200 µl fresh *E. coli* (recipient) o/n culture to each tube. Incubate 20 min at 37 °C.
- 4.** Pellet cells and resuspend in 1 ml LB + 20 mM sodium citrate. Incubate 1 hour at 37°C, no shaking.
- 5.** Pellet and resuspend in 1ml TM (10mM Tris HCL pH 7.5, 10 mM MgCl₂).
- 6.** Pellet and resuspend in 50 µl TM. Plate on selective media (antibiotic marker is being transferred by the P1 lysate). If the donor strain is *lac*⁻ and the recipient strain is *lac*⁺ plate onto MacConkey + lactose + antibiotic plates and select only red colonies.
- 7.** Verify strain by PCR.